

www.texwipe.com


Tex Mop<sup>™</sup> Mop Handles

H. SIGIL OH HILLING



# Cleanroom Mops TexMop<sup>™</sup> Handles

# TECHNICAL DATA SHEET



# **Products**

Number	Description	Sterile	Packaging	Weight
Works wit	h AlphaMop ™, Mini AlphaMop ™/Isolator Cleaning Tool, ClipperMop ™			
* TX7117	16" - 23" (41 cm - 58 cm), Telescoping, Stainless Steel		1 handle/case	2 lbs. (0.9 kg)
TX7120	12" (30 cm), Fiberglass, White		1 handle/case	2 lbs. (0.9 kg)
TX7121	18" (46 cm), Fiberglass, White		1 handle/case	2 lbs. (0.9 kg)
TX7122	29" - 53" (74 cm - 1.35 m), Telescoping, Fiberglass, White		1 handle/case	2 lbs. (0.9 kg)
TX7123	60" (1.5 m), Fiberglass, White		1 handle/case	2 lbs. (0.9 kg)
* TX7181	4' - 16' (1.2 m - 4.8 m), 4 piece Extension Handle, Telescoping, Aluminum, Blue		1 handle/case	5 lbs. (2.3 kg)
* TX7182	12" - 30" (30 cm - 76 cm), Telescoping, Aluminum, Blue		1 handle/case	2 lbs. (0.9 kg)
* TX7183	55" - 92" (1.4 m - 2.3 m), Telescoping, Aluminum, Blue		1 handle/case	2 lbs. (0.9 kg)
TX10217	TexConnectTM Adapter, Aluminum		1 adapter/case	1 lb. (0.5 kg)
Works wit	h BetaMop ™ II			
TX7125	60" (1.5 m), Fiberglass, White		1 handle & 1 clamp/case	2 lbs. (0.9 kg)
TX7092	60" (1.5 m), Stainless Steel		1 handle & 1 clamp/case	4 lbs. (1.8 kg)

\* Use with TX10217 TexConnect Adapter







## TECHNICAL DATA SHEET

## **Description**

Texwipe offers a complete line of mop handles that can be used interchangeably with our Cleanroom mops including AlphaMop<sup>™</sup>, BetaMop<sup>™</sup>II, and ClipperMop<sup>™</sup>.

# **Features & Benefits**

- Mop handle lengths range from 12 inches (30 cm) to 16 feet (4.8 m) in order to clean confined areas as well as hard-to-reach areas.
- The TexConnect<sup>™</sup> adapter allows for interchangeable use with most screw-based mop heads.
- Stainless steel, aluminum and fiberglass options
- Good chemical resistance for compatibility with a variety of solutions.
- All handles and the adapter are autoclave safe (20 minutes @ 250°F/121°C)
- Lightweight, ergonomic, stable design

### **Applications**

- Cleaning ceilings, floors, walls , isolators, hoods and other confined areas
- Applying and removing solutions and disinfectants

# **Cleanroom Environment**

- ISO Class 2 7
- Class 1 100,000
- EU Grade A D

### Industries

Aerospace	Animal Laboratory	Biologics
Cleanroom Design/Build	Compounding Pharmacies	Data Storage
Facilities Maintenance	Industrial	Laboratory
Medical Device	Microelectronics	Pharmaceutical
Printing/Graphics	Semiconductor	USP <797> / USP <800>







# TECHNICAL DATA SHEET

#### Mops

Number	Description		
AlphaMop™			
TX7108	AlphaMop™ Cleanroom Mop with head assembly and handle		
TX7108	AlphaMop <sup>™</sup> - Autoclavable - Cleanroom mop with head assembly and handle		
TX7108H	AlphaMop™ replacement head assembly only		
TX7108HA	AlphaMop™ - Autoclavable - replacement head assembly only		
Mini AlphaM	op™/Isolator Cleaning Tool		
TX7101	Mini AlphaMop™ Isolator Cleaning Tool with head assembly and handle		
TX7104	Mini AlphaMop™ with head assembly and handle		
TX7105	Mini AlphaMop™ Isolator Cleaning Tool replacement head assembly only		
BetaMop™ /	11		
TX7106	BetaMop™ II with clamp, head refill and handle		
ClipperMop <sup>1</sup>	M		
TX7102	ClipperMop <sup>™</sup> - 7" x 4" (17.8 cm x 10 cm) with head assembly and handle		
TX7103	ClipperMop™ - 11" x 4" (27.9 cm x 10 cm) with head assembly and handle		
TX7112	ClipperMop <sup>™</sup> - 7" x 4" (17.8 cm x 10 cm) replacement head assembly only		
TX7113	ClipperMop™ - 11" x 4" (27.9 cm x 10 cm) replacement head assembly only		

### Mop Accessories by Texwipe



Texwipe's complete line of Cleanroom TexMop™ products are available at www.texwipe.com.

#### Texwipe holds ISO 9001 and ISO 14001 registrations.

All Texwipe products conform to GHS classification for labeling. Shipping classification based on weight of inner package.